






How to use DS-11+ (Nucleic acid measurement)

- ※The minimum required volume is 0.5ul. **1ul or more is recommended.**
- ※Because the optical path length is constant, the resultant density does not change due to the difference in volume.
- ※Measure immediately after putting the sample in the measuring position.
- ※Low concentration samples take time to measure, but it is not a malfunction.

1. **Clean** the measuring position. **Turn on** the DS-11+.
2. From the home menu **tap** measurement mode. (**dsDNA , RNA , ssDNA**)
3. **Put a blank** solution and **tap 'Blank'**. **Wipe** blank solution.
(**Wipe both the top and bottom** of the arm.)
4. (If necessary input sample name.) **Put the sample** and **tap 'Measure'**.
Waveform and density results are displayed. **Wipe** the sample.
5. Repeat sample measurement.
Swipe the measurement screen to the left to see a list of results.
Swipe to the left again and the graph will be displayed.
Tap on the sample name to display the graph.
Swipe to the right to return.
6. How to save data to USB memory.
Select sample name in result table or graph.
Select '**Export Selected Samples**' or '**Screen Capture**' from the **button at the upper right of the screen.**
'**Export Selected Samples**' stores the csv file of the **result table.**
'**Screen Capture**' stores the jpeg file of the **graph.**
7. When all measurement is completed, select '**Exit**' from the **button at the upper right of the screen.**
Clean the measuring position.
8. **Confirm that write operation is not in progress** and disconnect USB.
9. Write the record.



Measurement screen (RUN)



Result table (REPORT)

Graph


How to use DS-11+ (Spectrum measurement)

- ※The minimum required volume is 0.5ul. **1ul or more is recommended.**
- ※Because the optical path length is constant, the resultant density does not change due to the difference in volume.
- ※Measure immediately after putting the sample in the measuring position.
- ※Low concentration samples take time to measure, but it is not a malfunction.

1. **Clean** the measuring position. **Turn on** the DS-11+.
2. From the home menu **tap** the **'UV-VIS'**.
3. Select the **cuvette type**. When you want to use by drop, select 'microvolume'.
4.  Select **'Add new'** from this icon. **Enter** the measurement **wavelength** in the lower row.
Enter the **'List name'** and tap the OK.
5. **Put a blank** solution and **tap 'Blank'**. **Wipe** blank solution.
(**Wipe both the top and bottom** of the arm.)
6. (If necessary input sample name.) **Put the sample** and **tap 'Measure'**.
Waveform and density results are displayed. **Wipe** the sample.
7. Repeat sample measurement.
Swipe the measurement screen to the left to see a list of results.
Swipe to the left again and the graph will be displayed.
Tap on the sample name to display the graph. Swipe to the right to return.
8. How to save data to USB memory.
Select sample name in result table or graph.
Select **'Export Selected Samples'** or **'Screen Capture'** from the **button at the upper right of the screen.**
'Export Selected Samples' stores the csv file of the **result table.**
'Screen Capture' stores the jpeg file of the **graph.**
9. When all measurement is completed, select **'Exit'** from the **button at the upper right of the screen.**
Clean the measuring position.
10. **Confirm that write operation is not in progress** and disconnect USB.
11. Write the record.



When AutoRun is checked, measurement is started when the arm is lowered.



Measurement screen (RUN)





Result table (REPORT)

Graph


How to use DS-11+ (Protein measurement)

- ※The minimum required volume is 0.5ul. **1ul or more is recommended.**
- ※Because the optical path length is constant, the resultant density does not change due to the difference in volume.
- ※Measure immediately after putting the sample in the measuring position.
- ※Low concentration samples take time to measure, but it is not a malfunction.


1. **Clean** the measuring position. **Turn on** the DS-11+.
2. From the home menu **tap** measurement mode. (**Protein A280** , **Labeled Proteins** , **Peptides**)
3. Select **sample type**.
4. **Put a blank** solution and **tap 'Blank'**. **Wipe** blank solution.
(**Wipe both the top and bottom** of the arm.)
5. (If necessary input sample name.) **Put the sample** and **tap 'Measure'**.
Waveform and density results are displayed. **Wipe** the sample.
6. Repeat sample measurement.
Swipe the measurement screen to the left to see a list of results.
Swipe to the left again and the graph will be displayed.
Tap on the sample name to display the graph.
Swipe to the right to return.
7. How to save data to USB memory.
Select sample name in result table or graph.
Select **'Export Selected Samples'** or **'Screen Capture'** from the **button at the upper right of the screen**.
'Export Selected Samples' stores the csv file of the **result table**.
'Screen Capture' stores the jpeg file of the **graph**.
8. When all measurement is completed, select **'Exit'** from the **button at the upper right of the screen**.
Clean the measuring position.
9. **Confirm that write operation is not in progress** and disconnect USB.
10. Write the record.



When AutoRun is checked, measurement is started when the arm is lowered.



Measurement screen (RUN)



Result table (REPORT)

Graph

How to use DS-11+ (Protein measurement by colorimetric method)


※The minimum required volume is 0.5ul. 1ul or more is recommended.

※Because the optical path length is constant, the resultant density does not change due to the difference in volume.

※Measure immediately after putting the sample in the measuring position.

※Low concentration samples take time to measure, but it is not a malfunction.

1. **Clean** the measuring position. **Turn on** the DS-11+.
2. From the home menu **tap** the 'Colorimetrics'.
3. Select **cuvette type** and **assay type**. When you want to use by drop, select 'microvolume'.
4. Displayed the 'STD Curve' screen when **swipe** to the left. (This is standard sample measurement mode)
5. Select 'New Curve' from the pulldown.
Tap the 'Std. Concentrations' and **enter standard concentrations**.
6. **Put a blank solution** and **tap 'Blank'**. **Wipe** blank solution.
(**Wipe both the top and bottom** of the arm.)
7. Put a standard solution. **Tap the standard concentration position** in the table.
(The row of ABS 2 represents the second measurement.) Tap the 'Measure'. **Wipe** the standard.
8. Repeat the standard measurement and **measure all standard solutions**. Confirm that the standard curve is correct.
9. Displayed the 'Run' screen when **swipe** to the right. (This is unknown sample measurement mode)
10. (If necessary input sample name.) **Put the sample** and **tap 'Measure'**.
Waveform and density results are displayed. **Wipe** the sample.
11. Repeat sample measurement.
Swipe the measurement screen to the left twice to see a list of results.
Swipe to the left again and the graph will be displayed. Tap on the sample name to display the graph.
Swipe to the right to return.
12. How to save data to USB memory.
Select sample name in result table or graph.
Select 'Export Selected Samples' or 'Screen Capture' from the **button at the upper right of the screen**.
'Export Selected Samples' stores the csv file of the **result table**. 'Screen Capture' stores the jpeg file of the **graph**.
13. When all measurement is completed, select 'Exit' from the **button at the upper right of the screen**.
Clean the measuring position.
14. **Confirm that write operation is not in progress** and disconnect USB.
15. Write the record.



How to use DS-11+ (Standard curve method)


※The minimum required volume is 0.5ul. 1ul or more is recommended.

※Because the optical path length is constant, the resultant density does not change due to the difference in volume.

※Measure immediately after putting the sample in the measuring position.


※Low concentration samples take time to measure, but it is not a malfunction.

1. **Clean** the measuring position. **Turn on** the DS-11+.
2. From the home menu **tap** the '**STD Curve Method**'.
3. Select the **cuvette type**. When you want to use by drop, select 'microvolume'.
4.  Select '**Add new**' from this icon.
Enter the measurement wavelength and baseline wavelength in 'Analysis nm' and 'Baseline nm' field.
Enter the **measurement range** in the 'min nm' and 'max nm' fields. Enter the '**Method Name**' and tap the OK.
5. Displayed the '**STD Curve**' screen when swipe to the left. (This is standard sample measurement mode)
Confirm that the method name is displayed on the screen.
Tap the '**STD. Concentrations**'. Enter the **concentration of the standard**.
6. **Put** a **blank** solution and **tap 'Blank'**. **Wipe** blank solution.(**Wipe** both the **top and bottom** of the arm.)
7. Put a standard solution. **Tap the standard concentration position** in the table.
(The row of ABS 2 represents the second measurement.) Tap the '**Measure**'. **Wipe** the standard.
8. Repeat the standard measurement and **measure all standard solutions**. Confirm that the standard curve is correct.
9. Displayed the '**Run**' screen when swipe to the right. (This is unknown sample measurement mode)
10. (If necessary input sample name.) **Put** the **sample** and **tap 'Measure'**.
Waveform and density results are displayed. **Wipe** the sample.
11. Repeat sample measurement.
Swipe the measurement screen to the left twice to see a list of results.
Swipe to the left again and the graph will be displayed. Tap on the sample name to display the graph.
Swipe to the right to return.
12. How to save data to USB memory.
Select sample name in result table or graph.
Select '**Export Selected Samples**' or '**Screen Capture**' from the **button at the upper right of the screen**.
'**Export Selected Samples**' stores the csv file of the **result table**. '**Screen Capture**' stores the jpeg file of the **graph**.
13. When all measurement is completed, select '**Exit**' from the **button at the upper right of the screen**.
Clean the measuring position.
14. **Confirm that write operation is not in progress** and disconnect USB.
15. Write the record.



How to use DS-11+ (How to save the data stored in the device)

1. **Turn on** the DS-11+.
2. From the home menu **tap** the 'Data'.
3. How to search data. **Select the measurement mode** in the Application field.
Tap the '**Search**' to display the corresponding data.
4. **Select** what you need from the list and tap the '**Add to Report**'.
5. Displayed result table or graph when **swipe** to left.
Tap the sample name on the graph screen to display the graph.
6. How to save data to USB memory.
Select sample name in result table or graph.
Select '**Export Selected Samples**' or '**Screen Capture**' from the **button at the upper right of the screen.**
'Export Selected Samples' stores the csv file of the **result table.**
'Screen Capture' stores the jpeg file of the **graph.**
7. **Confirm that write operation is not in progress** and disconnect USB.
8. Select '**Exit**' from the **button at the upper right of the screen.**




How to use DS-11+ (Kinetic)

※The minimum required volume is about 1 ml with standard cuvette.

※Because the optical path length is constant, the resultant density does not change due to the difference in volume.

1. **Clean** the measuring position. **Turn on** the DS-11+.
2. From the home menu **tap** the 'Kinetics'.
3. Select the **cuvette type**.
4.  Select 'New method' from this icon.
Enter the measurement wavelength and baseline wavelength in 'Analysis nm' and 'Baseline nm' field.
Enter the **measurement range** in the 'min nm' and 'max nm' fields.
Enter the **temperature** in the 'Set Heater Temp' field. (Uncheck when temperature control is not necessary.)
When using a quartz cuvette, check the 'Use Quartz Cuvette'.
Enter delay time from start in 'Delay before 1st Measurement' field.
Enter **interval time** in 'Measurement Interval' field.
Enter **measurement time** in 'Duration of Stage' field. Enter the 'Method Name' and tap the confirm button.
5. To control the temperature, turn on the heater control. Wait until the temperature stabilizes.
6. **Put** a blank solution and **tap 'Blank'**. **Wipe** blank solution. (Wipe both the **top and bottom** of the arm.)
7. (If necessary input sample name.) **Put** the **sample** and **tap 'Measure'**.
(A graph of measurement wavelength is displayed during measurement.)
8. (After measurement) A list of results is displayed on the report screen.
The spectrum of each point is displayed on the graph screen. A kinetic graph is displayed on the trend screen.
9. How to save data to USB memory.
Select sample name in result table or graph.
Select 'Export Selected Samples' or 'Screen Capture' from the **button at the upper right of the screen.**
'Export Selected Samples' stores the csv file of the **result table**. 'Screen Capture' stores the jpeg file of the **graph**.
10. When all measurement is completed, select 'Exit' from the **button at the upper right of the screen.**
Clean the measuring position.
11. **Confirm that write operation is not in progress** and disconnect USB.
12. Write the record.

Caution!
The display is stopped during the delay cycle.



Main switch

Insert cuvette

Kinetics
Mode selection

Method name: Analysis nm 1, Analysis nm 2
 Enter New Method Name: Min nm: 7300, Range: 700, Baseline nm: 700
 Set Heater Temp: Temperature control and temperature
 Use Quartz Cuvette: Quartz Cuvette
 Stage 1 Stage 2 Stage 3
 Delay before 1st Measurement: Delay time: 00:00:10
 Measurement Interval: Interval time: 00:00:15
 Duration of Stage (including delay): Measurement time: 00:01:00
 Note: Durations are rounded down to the nearest non-multiple of intervals, so 31 seconds at 10 second intervals is 30 seconds.

Blank Measure Mode
 Sample name:
 Savefile: 10.mpg
 Temperature control:

Export, Screen capture, Exit
In this button

全選択
swipe
Select sample
swipe
swipe

Measurement screen (RUN) Result table (REPORT) Graph Trend