Okayama University Medical Research Updates (OU-MRU) 2016.03 Vol.22

Source: Okayama University (JAPAN), Public Relations and Information Strategy

For immediate release: 22 March 2016

Okayama University research: Medical supportive device for hemodialysis catheter puncture

(Okayama, 22 March 2016) Medical doctors at Okayama University Medical School and Shigei Medical Research Hospital in Okayama City, in collaboration with K.Techno Inc. develop a medical supportive device for hemodialysis catheter puncture.

Hemodialysis is used to purify the blood of people with kidney failure with a growing need for dialysis using a catheter in countries with aging populations, such as Japan. The precise catheter placement for secure and efficient dialysis, requires specialized knowledge and skills of the doctor.

The implication of improper catheter placement include: the appropriate and correct hemodialysis cannot be carried out; the possibility of medical accidents;, increase the exposure of patients to radiation through the use of X-ray equipment during replacement; and wastefulness in healthcare costs due to unnecessary disposal of expensive dialysis catheters.

Now, Dr Toshiaki Ohara at the Okayama University Medical School and clinical staff at the Shigei Medical Research Hospital in Okayama City, in collaboration with K.Techno Inc., have developed a "medical supportive device for hemodialysis catheter puncture" (hereafter 'device').

Clinical tests conducted at Shigei Medical Research Hospital showed no safety problems with excellent placement results of hemodialysis catheters.


The Okayama team will publish their results in an international journal and are preparing to commercialize this invention as well as looking for other non- hemodialysis catheter related applications of this technology.

Details of the process for using the 'device'

Before inserting the catheter the patient lies horizontally on an X-ray table followed by disinfection and sheets to define the operation area. Next, the 'device' is laid over the patient and a pen used to mark the path for the catheter. Then, a catheter of the appropriate length is selected and the internal jugular vein is punctured.

Then, in a state in which the wire and the dilator have been inserted, the device is used to make a line along the wire again. Then, actually superimpose the catheter to check that there is no deviation.


Then the exit section is set and the catheter is guided in a direction from the outlet to the neck, and the inner cylinder of the dilator is removed. The catheter is inserted into the dilator and the process is completed when the dilator sheath is split and removed. This process enables the accurate and reproducible placement of catheters without deviation.



Medical supportive device for hemodialysis catheter puncture




Patient and medical supportive device for hemodialysis catheter puncture



X-ray images before (left) and after placement of catheter

Correspondence to

Assistant Professor Toshiaki Ohara, M.D., Ph.D. Department of Pathology & Experimental Medicine, Graduate School of Medicine, Dentistry and Pharmaceutical Sciences, Okayama University, Shikata-cho 2-5-1, Okayama city, Okayama 700-8558, Japan E-mail: t_ohara@cc.okayama-u.ac.jp



Assistant Professor Toshiaki Ohara

Further information

Okayama University 1-1-1 Tsushima-naka , Kita-ku , Okayama 700-8530, Japan Public Relations and Information Strategy E-mail: www-adm@adm.okayama-u.ac.jp Website: http://www.okayama-u.ac.jp/index e.html Okayama Univ. e-Bulletin: http://www.okayama-u.ac.jp/user/kouhou/ebulletin/ Okayama Univ. e-Bulletin: http://www.okayama-u.ac.jp/user/kouhou/ebulletin/ Okayama Univ. e-Bulletin (PDF Issues): http://www.okayamau.ac.jp/en/tp/cooperation/ebulletin.html About Okayama University (You Tube): https://www.youtube.com/watch?v=iDL1coqPRYI Okayama University Image Movie (You Tube): https://www.youtube.com/watch?v= WnbJVk2eIA https://www.youtube.com/watch?v=KU3hOIXS5kk Okayama University Medical Research Updates (OU-MRU) Vol.1 : <u>Innovative non-invasive 'liquid biopsy' method to capture circulating tumor cells</u> from blood samples for genetic testing

- Vol.2 : Ensuring a cool recovery from cardiac arrest
- Vol.3 : Organ regeneration research leaps forward
- Vol.4 : Cardiac mechanosensitive integrator
- Vol.5 : <u>Cell injections get to the heart of congenital defects</u>
- Vol.6 : Fourth key molecule identified in bone development
- Vol.7 : Anticancer virus solution provides an alternative to surgery
- Vol.8 : Light-responsive dye stimulates sight in genetically blind patients
- Vol.9 : Diabetes drug helps towards immunity against cancer
- Vol.10 : Enzyme-inhibitors treat drug-resistant epilepsy
- Vol.11 : Compound-protein combination shows promise for arthritis treatment
- Vol.12 : Molecular features of the circadian clock system in fruit flies
- Vol.13 : Peptide directs artificial tissue growth
- Vol.14 : Simplified boron compound may treat brain tumours
- Vol.15 : Metamaterial absorbers for infrared inspection technologies
- Vol.16 : Epigenetics research traces how crickets restore lost limbs
- Vol.17 : Cell research shows pathway for suppressing hepatitis B virus
- Vol.18 : Therapeutic protein targets liver disease
- Vol.19 : Study links signalling protein to osteoarthritis
- Vol.20 : Lack of enzyme promotes fatty liver disease in thin patients
- Vol.21 : Combined gene transduction and light therapy targets gastric cancer

About Okayama University

Okayama University is one of the largest comprehensive universities in Japan with roots going back to the Medical Training Place sponsored by the Lord of Okayama and established in 1870. Now with 1,300 faculty and 14,000 students, the University offers courses in specialties ranging from medicine and pharmacy to humanities and physical sciences.

Okayama University is located in the heart of Japan approximately 3 hours west of Tokyo by Shinkansen.

